

Kinreen Technology (Shenzhen) Co.,Ltd

Company Profile and Product Catalogue

FDA Registered LED Light Manufacturer CE, FCC,RoHS, SGS Certification

OEM & ODM & Drop-shipping & Fast Shipping

Website:<https://www.kinreen.com>

— ABOUT US

Kinreen Tech. (Shenzhen)Co., Ltd. is a professional Alibaba and SGS verified factory manufacturing Led therapy lights .We have more than 17 years experience with well-equipped R & D capabilities and excellent management to guarantee quality product .

Kinreen's values :

Company mission : Develop products for happiness life .

Company prospect : Make life more enjoyable .

Company core values : Honesty and trustworthiness ; Fairness ; Innovation;

Main market : North America , Europe , Oceania ;

All of our products are CE , ROHS , FCC approved for business all over the world . We are doing business with up to 897 clients from all over the world, which are from US, Canada, Australia, New Zealand , UK, France , Germany, Italy, Austria, Thailand, Singapore, Brazil, South Africa, etc.

Kinreen pay much attention to the needs of the customer , rely on excellent product quality , reasonable price and fast delivery , won the recognition and trust of customers step by step, which can be your reliable supplier, and will always be.



— QUALITY CONTROL SYSTEM —



Production Inspection



Aging Testing



Customized Your Brand



Inspection Before Shipment



Testing Data Before Shipment



Packing

CERTIFICATE

Certificate • Сертификат • 證書 • Certificat • 증명서 • 證明書

CERTIFICATE OF CONFORMITY

Registration No. WST19N030058
Technical Construction File No. WST19N030058-1ER

Certificate Holder: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Manufacturer: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Product: LED Therapy Light
Model (S): QY-S-B04, QY-S-C74, QY-S-E24, QY-S-E26, QY-P-D04, QY-P-D06, QY-P-A0-C, QY-P-A0-C

Trade Mark: **Kinreen**

Rating: 100-240VAC, 50/60Hz, 30W

Directives: EMC, LVD, RoHS

Standard(s): GB 4793.1-2015, GB 4793.2-2015

RoHS

**Fiscal Year 2020
CERTIFICATION OF REGISTRATION**

This certifies that:

KINREEN TECHNOLOGY(SHENZHEN) CO., LTD.
207-208 Chuangke Building 72-6 Huangnan South Road Xintian Community Guanta Street, Longhua District Shenzhen, Guangdong, 518110, CHINA

has completed the FDA Establishment Registration (as manufacturer, contract manufacturer) and Device Listing with the US Food & Drug Administration, through

U.S. Agent for FDA Communications: **SUNGO TECHNICAL SERVICE INC.**
6950 W EASTWOOD AVE APT 201, CHICAGO, ILLINOIS 60630, USA
Telephone: +1-855-937-7779 / E-mail: sungo.group@yahoo.com

Owner/Operator Number: 10061531

Device Listing#:

Listing No	Code	Device Name
D-1-5	ELY	Lamp, infrared, therapeutic heating (LED Therapy Light LED)

SUNGO Technical Service Inc. will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above, unless said registration is terminated after issuance of this certificate. SUNGO Technical Service Inc. makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not constitute endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. SUNGO Technical Service Inc. assumes no liability to any person or entity in connection with the foregoing.

Paragraph 31 CFR 807.18: "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding." The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. SUNGO Technical Service Inc. is not affiliated with the U.S. Food and Drug Administration.

Sungo FDA

Executive Director
Issued: Oct. 28 2019
Cert. No. 200815362518
Expiration Date: Dec. 31 2020

SUNGO CHINA OFFICE Tel: 021-68828052 Email: Shago2008@126.com Website: www.sungoiprimal.com
Add: 12th Floor, No.1500 Century Avenue, Shanghai 200122, P.R.China

VERIFICATION OF CONFORMITY

According to FCC Part 15B

Registration No. WST19N030053E
Technical Construction File No. WST19N030053-1ER

Certificate's Holder: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Manufacturer: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Product: LED Therapy Light
Model (S): QY-S-B04, QY-S-C74, QY-S-E24, QY-S-E26, QY-S-E28, QY-P-D04, QY-P-D06, QY-P-A0-C, QY-P-A0-C

Trade Mark: **Kinreen**

Rating: 100-240VAC, 50/60Hz, 30W

Standard(s): FCC Part 15B Subpart B
Procedure ANSI C63.14 2014

Compliance with Part 15B of FCC Rules
This Certificate with Part 15B of the FCC Rules. Operation is subject to the following two conditions: (1) The device may not cause harmful interference, and (2) The device must accept any interference received, including interference that may cause unintended operation.

Responsible Party: _____ Tested By: _____
Responsible Signature: _____ Issued By: _____
Date: _____ Date: Apr. 04, 2019

The Certificate of Verification shows that the tested sample technically complies with the FCC Part 15B. The certificate applies to the exact sample shown mentioned and should not imply an assessment of the entire.

CERTIFICATE OF CONFORMITY

Registration No. WST19N030052E
Technical Construction File No. WST19N030052-1ER

Certificate's Holder: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Manufacturer: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Product: LED Therapy Light
Model (S): QY-S-B04, QY-S-C74, QY-S-E24, QY-S-E26, QY-S-E28, QY-P-D04, QY-P-D06, QY-P-A0-C, QY-P-A0-C

Trade Mark: **Kinreen**

Rating: 100-240VAC, 50/60Hz, 30W

Standard(s): GB 4793.1-2015, GB 4793.2-2015

RoHS



VERIFICATION OF CONFORMITY

Registration No. WST19N030053E
Technical Construction File No. WST19N030053-1ER

Certificate Holder: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Manufacturer: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Product: LED Therapy Light
Model (S): QY-S-B04, QY-S-C74, QY-S-E24, QY-S-E26, QY-S-E28, QY-P-D04, QY-P-D06, QY-P-A0-C, QY-P-A0-C

Trade Mark: **Kinreen**

Rating: 100-240VAC, 50/60Hz, 30W

Standard(s): FCC Part 15B Subpart B
Procedure ANSI C63.14 2014

United States of America
Kinreen

Big The Label
Registered with the FDA

Manufacturer: Kinreen Technology (Shenzhen) Co., Ltd.
1558 Changhe Building 72-6 Huangnan South Road
Guanta Street Longhua District Shenzhen

Product: LED Therapy Light
Model (S): QY-S-B04, QY-S-C74, QY-S-E24, QY-S-E26, QY-S-E28, QY-P-D04, QY-P-D06, QY-P-A0-C, QY-P-A0-C

Trade Mark: **Kinreen**

Rating: 100-240VAC, 50/60Hz, 30W

Standard(s): GB 4793.1-2015, GB 4793.2-2015

外观设计专利证书

外观设计专利证书

专利号: 201930123456.7

申请人: 深圳金润科技有限公司

发明人: 张三, 李四

专利权人: 深圳金润科技有限公司

专利代理机构: 深圳金润知识产权代理有限公司

公告日期: 2019年12月31日

专利权期限: 十年



专利地址: 广东省深圳市龙岗区



protect your payment secure your order

100%
SECURITY GUARANTEED

Hot Sale Red Light Therapy Wraps Collection



RED LIGHT THERAPY PAD



STRETCHY BAND



ADAPTOR



POWER CORD



CONNECTOR WIRE

KR120 60W

Size: 400*200mm



KR200 100W

Size:492*238mm




KR75 37.5W

Size:1160*104mm

Newest Red Light Therapy Wraps

(built-in battery ,no need to charge all the time)



KR30 Wrap

KR72 Wrap

KR 64 Wrap

Large Size Wraps for Full Body , For Horses




KR180 Wrap one pair=2 pads

Size:800*160*6.3mm
Led pcs: 180*2 pcs 0.5W

540 chips of each pad
Wavelength: 660nm:850nm=1:2

All our wrap series is made of SBR waterproof material, which is easy to clean and most important for quality assurance for one year



KR350 WRAP

Size: 890*280mm
Led pcs :350 pcs 0.5W


1050 chips of each pad
Wavelength: 660nm:850nm=1:2

Application : for Pain Relief on Body, Knee, Ankle, Hands, Feet, Dogs, Horses

Don't fold the device over 90 degree as below photo shows otherwise it may damage the wrap . Please kindly noted that any product problems cause by incorrect use of our products are not covered by our warranty policy .



Please wrap the device to your body correctly as below photo guides



Plug the power bank , so you can walk around while working and doing exercise .



Red Light Therapy Cap for Hair Regrowth



The device is ideal for hair re-growth ;
memory improvement ; Headache pain
relief ; Parkinson ; Alzheimer's treatment

150 PCS LED lights with 450 LED Chips

Wavelength: 630nm:850nm:940nm



Timer: 15 Minutes . The device
will shut off automatically after working for
15 minutes ;

Pulse mode :10HZ and 40 HZ

Pulse which show much more effects than
constant lighting



It can be equipped with power bank(Not included). You can walk
around freely and get rid of the fixed power supply .

Handheld 3LEDs & 5 LEDs red light therapy torch

Application : small treatment joint pain relief , hands, knee, ankle (human and pets)



3 LEDs torch



5 LEDs torch



5 LEDs torch with the oral

Model	Led pcs	Power	Color	Material	Wavelength	Timer	Pulse
3 leds torch	3 PCS	9W	Silver	Stainless steel	630nm:660nm:850nm	5 Minutes	No
5 leds torch	5 PCS	15W	Champaign gold	Aluminum	460nm:630nm:660nm:850nm:940nm	3 minutes	10HZ
5 leds torch with oral	5 pcs	15W	Champaign gold	Aluminum	460nm:630nm:660nm:850nm:940nm	3 minutes	10HZ

Red Light Therapy Slipper & Shoes For The Foot



Slipper

88 LEDs of each



Shoe

108 LEDs of each



Boot

296 LEDs of each

2023 NEW Arrive Therapy Pads



KR120 (Pad and strap in a whole)

KR160 (Built in battery)

KR120 (Built in Battery)

12 LEDs red light therapy



PRODUCT DETAILS



3 LIGHTING MODES



All LED lighting

Red light

Infrared light mode
(Invisible infrared light)



BETTER BEAUTY, MORE HEALTH



15 LEDs red light therapy




PRODUCT DETAILS



Infrared light indicator
Switch of the product /
Adjustment button



Power display
TYPE C charging port



3 LIGHTING MODES


15 min Automatic timing for 15 minutes

All LED lighting

Red light

Infrared light mode
(Invisible infrared light)

Led Red Light Therapy Face Mask



240 LED beads in total

Four light sources can be switched



Red light: 630nm



Blue light: 460nm



Yellow light: 590nm



Infrared light: 850nm

(Invisible to the naked eye, it can be opened at the same time with other lights)

	<p>Wear to work No need to hold hands</p>
	<p>For entertainment Beauty while skin care</p>
	<p>Do as one wishes No need to close your eyes No need to plug in</p>

Led Red Light Therapy Face Mask + Neck Mask

Feature : Acne Treatment, Skin Rejuvenation, Wrinkle Remover, Skin Tightening



Four kinds of light wave depth repair

460nm Blue light	590nm Yellow light	630nm Red light	850nm Infrared light
Sensitive repair Balance grease	Brighten skin Improve darkening	Enhance elasticity Resist fine lines	Repair skin blood circulation


The diagram shows a cross-section of the skin with four vertical lines representing light waves of different colors and wavelengths. The blue line (460nm) is the shallowest, followed by yellow (590nm), red (630nm), and infrared (850nm) being the deepest. The infrared line is shown reaching the deepest layer of the skin.

Double repair of mask neck

Ensure stereo coverage of light source
No choice of face shape to fit the face

The image shows a person wearing the mask with a neck band. To the right are four icons with text: a water drop icon for "Enhanced absorption", a circular arrow icon for "Resist early aging", a wavy line icon for "Reduce wrinkles", and an upward arrow icon for "Repair skin".

Eye Massager



Facial Beauty Device




Product parameters



160 times / second high-frequency vibration massage

Relaxes pores Accelerates blood circulation to the face








360° free rotating guide head

Free guide head "zero gap" to fit the curve of the face




KR2400&KR3600 Led Red Light Therapy Panels



		
Model	KR2400W panel	KR3600W panel
Led Power	2800W	4200W
Led Chip	3W/5W*560pcs	3w/5w*840pcs
Wavelengths	5 wavelengths	
Red & Near Infrared	630:660:810:830:850nm	
Controller	Remote Control, Digital Display, 3 light Model, Timer (0-30) ,5 Dimmer , Pulse (10-587HZ)	
Irradiance	>190mw/cm2	>190mw/cm2
Light Angle	30 degree	
Dimensions	45.7"*16.5'*2.8'	67.7"*16.5"*2.8"
Net Weight	48lb	66lb
Lifespan	1000000 hrs	
EMF	0.0ut	
Best For	Face, Full of Body	
Warranty	1 year	

Folded 2/3/4/5 Therapy Panels





FACIAL BEAUTY



PROMOTES WOUND HEALING




MUSCULAR RELAXATION



WHOLE BODY USE

High Power Full Body Use , Face Beauty , Skin Rejuvenation Panels



KR100W



KR300W



KR300W with stand



KR1000W

Model	Power	Led pcs	Wavelength	Timer	Pulse	Function
KR100 W	100W	20 pcs 5W	660nm:850nm	10/20/30 Minutes	10HZ/40HZ	Pain relief /skin care
KR300W	300w	60PCS 5W	660nm:850nm	10/20/30 Minutes	10HZ/40HZ	Pain relief /skin care
KR300W with stand	300w	60PCS 5W	660nm:850nm	10/20/30 Minutes	10HZ/40HZ	Pain relief /skin care
KR1000W	1000w	200 PCS 5 W	660nm:850nm	10/20/30 Minutes	10HZ/40HZ	Pain relief /skin care

Different Size Red light Therapy Belt for Choice

Model	Size	LED qty	Power	Input AC with Adapter	Controller	Wavelength	Timing function	Pulse model	Power bank supply	Application
KR120	400*200mm	120	60W	YES	YES	660nm:850nm=1:2	20 Minutes	10HZ	YES	Waist, back, belly,pets
KR75	1160*104mm	75	37.5W	YES	YES	660nm:850nm=1:2	20 Minutes	10HZ	YES	Shoulder, waist,pets
KR200	492*238mm	200	100W	YES	YES	470nm:660nm:850nm	20 Minutes	10HZ	YES	Face ,waist,
KR350	890*280mm	350	175W	YES	YES	660nm:850nm=1:2	20 Minutes	10HZ	YES	Horses , pets, full body ,horses
KR180	880*160mm	180	90W	YES	YES	660nm:850nm=1:2	20 Minutes	10HZ	YES	Foots, waist ,waist
KR72		72	36W	Built-in battery		660nm:850nm=1:2	15 Minutes	10HZ		Knee, wrist, arms, leg
KR64	590*1650 mm	64	32W	Built-in battery		660nm:850nm=1:2	15 Minutes	10HZ		Neck, wrist
Kr30	550*7 0mm	30	15W	Built-in battery		660nm:850nm=1:2	15 Minutes	10HZ		Neck ,head

Red Light Therapy Main Function

Effectively remove wrinkles and fine lines



Good for hair growth

Reduce age spots and under-eye circles




Promote wound healing



Slimming




Relieves joint pain



Effective treatment for acne
Blackhead remove



HOW RED LIGHT THERAPY WORKS



As we age, we produce less collagen which results in the appearance of fine lines and wrinkles due to thin, dry skin.

Red Light Therapy penetrates deep into the dermis layer, stimulating cells that produce collagen.

New collagen moves to the surface of the skin, reducing fine lines and wrinkles.

The increase in collagen production smoothes the skin and reduces the appearance of fine lines and wrinkles.

Thanks for your time ! Please feel free to contact us for any questions .