

# ELECTROTHERAPY DEVICES CATALOGUE

---





# R-C101A

## (TENS + EMS + RUSS + IF)

- 2** Channel
- 4** IN 1 (TENS+EMS+RUSS+IF)
- 4** 50\*50mm pads (Japanese gel)
- 10** Body parts
- 90** Levels
- 90** mA at 1000 Ohm load
- 100** Programs
- 1050** mAh Rechargeable Li-ion battery




### SPECIFICATION

LCD Type	HTN
Backlight	white
Color	Black(optional)
Treatment frequency	2Hz-120Hz
Treatment pulse width	50uS-400uS
Treatment time	5-90min adjustable
Waveform	Symmetrical Bi-phase square-wave pulse
Dimensions	120.5x69.5x27mm
IP rating	IP22


### STIMULATION FORMS CONVENTIONAL



Burst



Frequency Modulated



Pulse Duration Modulated



Frequency And Pulse Duration Modulated

### PROGRAMS

#### TREATMENT AREA

#### + 10 body parts

Neck, Shoulder, Arm, Hand, Abdomen, Back, Hip, Leg, Foot, elbow, knee, ankle

#### TREATMENT PROGRAM

- + TENS: 30 programs** (Includes 10 customizable programs)
- + EMS: 20 programs** (Includes 7 customizable programs)
- + RUSS: 20 programs** (Includes 7 customizable programs)
- + IF: 30 programs** (Includes 10 customizable programs)

# TREATMENT AREA




Shoulder



Back



Waist



Abdomen



Foot



Leg



# R-C101C (4CH TENS+EMS)

## PRODUCT FEATURES

- 2** IN 1 (TENS+EMS)
- 4** Channel
- 8** 50\*50mm pads (Japanese gel)
- 12** Body parts
- 40** Levels
- 120** Max. 120mA (at 500 Ohm load)
- 1050** mA Rechargeable Li-ion battery



## SPECIFICATION


LCD Type	HTN
Backlight	white
Color	Black(optional)
Treatment frequency	2Hz- 120Hz
Treatment pulse width	50uS-400uS
Treatment time	30mins(5-90min adjustable)
Waveform	Symmetrical Bi-phase square-wave pulse
Dimensions	144x86x29.6mm
IP rating	IP22






# R-C101C (4CH TENS+EMS)


## STIMULATION FORMS CONVENTIONAL



Burst



Pulse Duration Modulated



Frequency Modulated



Frequency And Pulse Duration Modulated

## PROGRAMS

### TREATMENT AREA

#### + TENS: 12 body parts

Neck, Shoulder, Arm, Hand, Abdomen, Back, Hip, Leg, Foot, elbow, knee, ankle

#### + EMS: 9 body parts

Neck, Shoulder, Arm, Hand, Abdomen, Back, Hip, Leg, Foot

### TREATMENT PROGRAM

#### + TENS: 30 programs (Includes 12 customizable programs)

#### + EMS: 20 programs (Includes 9 customizable programs)

Optional Adjustable Program



# F100 (Foot Massager)

## PRODUCT FEATURES

**2** IN 1 (FOOT MASSAGER+BODY EMS)

**2** Channel

**2** 50\*50mm pads (Japanese gel)

**90** Levels

**90** Max. 90mA (at 1000 ohm)

**1050** mA Rechargeable Li-ion battery



## SPECIFICATION

LCD Type	HTN
Backlight	white
Color	Black(optional)
Treatment frequency	10Hz-36Hz
Treatment pulse width	250uS
Treatment time	10-90min adjustable
Waveform	Symmetrical Bi-phase square-wave pulse
Dimensions	367x361x80.5mm
IP rating	IP22



**F100 (FOOT MASSAGER)**

**INTENDED USE**



Foot



Leg Massages



Muscle Exercises

**POWER SUPPLY**

- ⚡ FOOT MASSAGER: D.C. 3.7V Li-ion Battery.
- ⚡ Remote Controller: D.C. 3.0V AAA Battery / 2pcs



**ACCESSORIES**

- ⦿ Foot massager-----1PC
- ⦿ Remote controller-----1PC
- ⦿ Electrode wire -----1PC
- ⦿ Electrode pad (50mm x 50mm)----2PCS
- ⦿ USB cable -----1PC
- ⦿ User manual -----1PC
- ⦿ Battery -----2PCS




# Mini TENS (with 4 designs Electrode)



## PRODUCT FEATURES

- ⊕ One treatment program and 4 treatment stages
- ⊕ Voice prompt function Timer
- ⊕ Function Pulse rate: 1Hz-100Hz
- ⊕ Pulse Width 100uS
- ⊕ Rechargeable Li-ion battery
- ⊕ Battery Includes USB Charging Cable



## MINI-TENS-COMBO Treatment Program list

Mode	Treatment Phase	Pulse Rate(Hz)	Pulse width(uS)	Phase time(Min)	Treatment time(Min)
TENS	1	1-100	100	3	20
	2	50	100	2	
	3	1	100	3	
	4	50/100	100	2	



# Mini TENS (with 4 designs Electrode)



Common Use



Hip Use



Abdomen Use



Back Use



# R-C101I (TENS)

- 2** Channel
- 4** 40\*40mm pads (Japanese gel)
- 5** Adjustable program  
(Burst , Normal, Modulation, Strength Duration 1, Strength Duration 2)
- 7** Preset programs
- 9** Volt Battery
- 60** Sets of operation records
- 100** Max.100mA (at 500 Ohm load)
- 999** Total treatment time up to 999 hours can be stored.



## SPECIFICATION

LCD Type	TN
Treatment frequency	2Hz-150Hz
Treatment pulse width	50uS-300uS
Treatment time	1-60mins and "C"
Dimension	101x61x24.5 mm (L x W x T)(without Belt clip)
Weight	About 150g
Carton	470*405*426mm, 32pcs/carton



## PROGRAMS

### Treatment Mode

+ TENS

### Treatment program

+ TENS: 12 programs  
(5 adjustable programs + 7 preset programs)

### Optional Adjustable Program

+ Width/Rate/Time



# M102A

- 2** Channel
- 3** IN 1 (TENS+EMS+MASSAGE)
- 4** 50\*50mm pads (Japanese gel)
- 12** Body parts
- 73** Programs
- 40** Levels
- 120** Max. 120mA (at 500 Ohm load)
- 500** mAh Rechargeable Li-ion battery



## SPECIFICATION

LCD Type	HTN
Treatment frequency	2Hz-120Hz
Treatment pulse width	50uS-300uS
Treatment time	5-90mins
Dimension	120x60x18.6mm (L x W x T)(without Belt clip)
Weight	About 90g
Carton	490*350*350mm, 50pcs/carton



## PROGRAMS

### Treatment Mode

+ TENS/EMS/MASSAGE

### Optional Adjustable Program

+ Hz/Us/Time

### Treatment program

+ TENS: 11 programs  
(adjustable 1 program)

+ EMS: 10 programs  
(adjustable 1 program)

+ MASSAGE  
(52 massage manipulations)





# R-C4A (Smart, fashionable)

- 2** Channel
- 3** IN 1 (TENS+EMS+MASS)
- 4** 50\*50mm pads (Japanese gel)
- 10** Body parts
- 40** Levels
- 60** Programs
- 120** Max.120mA (at 500 Ohm load)
- 500** mAh Rechargeable Li-ion battery

### BODY PARTS

TENS: Neck, Shoulder, Arm, Hand, Back, Abd.hip, Leg, Foot, Joint.

EMS: Neck, Shoulder, Arm, Hand, Back, Abd. Hip, Leg, Foot.

MASSAGE: Knead, Rub, Tap



## SPECIFICATION

LCD Type	HTN
Backlight	white
Dimension	109*54.5*23mm (L x W x T)(without Belt clip)
Weight	about 82g (without battery)
Carton	490*350*350mm, 50pcs/carton 13kgs/carton
Treatment channel	2
Treatment strength	Max.120mA (500 Ohm load)
Treatment intensity	0-40 levels,adjustable
Treatment frequency	2Hz-120Hz
Treatment pulse width	50uS-300uS
Treatment time	5-90mins



## PROGRAMS

### Treatment Mode

+ TENS/EMS/MASSAGE

### Optional Adjustable Program

+ Hz/Us/Time

### Treatment program

+ TENS: 30 programs (adjustable 10 programs)

+ EMS: 27 programs (adjustable 9 programs)

+ MASSAGE (3 massage manipulations)



## Belt Accessories



Care-Joint Belt  
for joint pain relief

AC3001



EMS ABS Fitness Belt  
for abdomen muscle training

AC1001



Care-Back Belt  
for back pain relief

AC2001

## Other hot sale Accessories



AC4001  
Gloves



AC5001  
Care-foot Socks



AC3002/AC3003  
Care-Knee(Elbow)  
Support



Clamp-A  
Suitable for R-C4  
series devices



Clamp-B  
Suitable for R-C4  
series devices



Pads Holder  
Suitable for 40\*40mm,  
50\*50mm,40\*80mm,  
50\*100mm Electrode pads

## Accessories: Lead wire



R-W101  
Lead Wire with  
Snap Connector



R-W102  
2-pin lead Wire  
(2.0mm)



R-W103  
4-pin lead Wire  
(2.0mm)



R-W104  
Lead Wire with  
2-Snap



R-W105  
2-pin lead Wire  
(2.0mm)




R-W106  
2-pin lead Wire  
(2.0mm)

## Accessories: Self-adhesive Electrode pads




40\*40mm/50\*50mm



40\*40/50\*50mm  
with Snap



40\*80mm/50\*100mm  
50\*130mm



40\*80mm EMS



32mm/40mm/50mm



40\*40/50\*50mm  
pads