



PRODUCT CATALOG

XIAMEN WEIYOU INTELLIGENT TECHNOLOGY CO.,LTD

Add:No.6 Xianghong Road, Xiang 'an District, Xiamen City, Fujian Province, China

Tel:(86)-592-6251545 Fax:(86)-592-6252548

Web:www.weiuit.com www.sunlionrecovery.com

E-mail:sales@weiuit.com



XIAMEN WEIYOU INTELLIGENT TECHNOLOGY CO.,LTD



MISSION VISION

Maintain life baseline and optimize life character. Let everyone can own a suitable recovery device.

CONTENTS

Company Profile	01-02
Warm Cold Compression Therapy Machine	03-04
Air Pressure Therapy System	05-28
Hand Rehabilitation Training Device	29-30
<hr/>	
Anti-Bedsores Inflatable Mattress With Pump	31-32
Medium Frequency Massager	33-34
Migraine Headache Relief	35-36
<hr/>	
Medical Parallel Transport Bed	37-38
Electric Medical Bed/Nursing Bed	39-40
Roller Massager	41-42
Product Display	43-46

COMPANY PROFILE

Xiamen Weiyu Intelligent Technology Co., Ltd. was established in early 2010, focusing on providing customers with comprehensive services of rehabilitation treatment, medical care equipment and intelligent consumer electronics. The team has rich industry experience in rehabilitation health care and intelligent technology, as well as technical R&D and innovation capabilities. Taking "high standard, high quality, high efficiency, high safety and high performance" as the standard, it strives to explore the development direction of health care life, and is determined to build a cross-industry ecological circle in health care field, so as to provide an overall solution of "maintaining life baseline and optimizing life character" for global health care life.

The company's main products include medical parallel transport bed series, air pressure therapy system series, anti-bedsore inflatable mattress with pump series, hand rehabilitation training device series, medium-frequency massager series, etc., and were certified by MFI (Made for iPhone) in the United States in 2012. With excellent products, reasonable prices, novel designs, environmentally-friendly materials and sincere service, our products have won the trust of customers at home and abroad, and our products are exported to Europe, America, the Middle East, Southeast Asia and other countries and regions. Become one of the most influential and competitive professional brands in the same industry.



Machine shop



WEIYOU INTELLIGENT TECHNOLOGY

Maintain life baseline and optimize life character.
Let everyone can own a suitable recovery device.



• Warm Cold Compression Therapy Machine

• Main purpose

Use cold compress, hot compress, and pressurization in one; relieve the pain of patients and improve blood circulation. For patients with long-term exercise, labor, various trauma, bone injury, fever, pain and swelling, etc., it can continuously improve the blood of the human body circulation and vasodilation.

It has the advantages of safety, simplicity, and high efficiency. It is an effective treatment method to treat chronic pain, promote blood circulation and remove blood stasis, relieve local muscle tension, and enhance immunity.

It can reduce swelling and pain in the early postoperative period, improve joint mobility, reduce complications, and bring unexpected therapeutic effects to patients.



• Warm Cold Compression Therapy Machine

▶ Machine Parameters

Voltage: DC12V	Time: Optional treating time
Power consumption: 400W	5~99 minutes
Weight: 10kgs	Pressure Range: 30-120mmHg
Size: 50*28*28.5cm	Water temperature range: 5°C~40°C



• Shoulder

• Waist

For different body parts



• Knee



• Foot + Calf



• Arm



Air Pressure Therapy System series

The main purpose

- The device is intended for use by healthcare professionals to administer treatment to patients or for patients to use at home by themselves. It is used to treat symptoms such as primary lymphedema, edema following trauma and sports injuries, post-immobilization edema, venous insufficiency, and lymphedema.

Applications in the medical field

- Prevention of deep vein thrombosis (DVT), venous insufficiency
- Various primary, secondary and mixed edema such as lymphedema and postpartum edema
- Fracture, femoral head necrosis, soft tissue injury and other bone injuries
- Peripheral neuritis and diabetic foot caused by diabetes
- Rehabilitation patients after surgery
- Paralyzed
- Patients with varicose veins
- Long-time standing, sedentary or bedridden patients

Applications in sports field

- Relax and relieve fatigue after exercise
- Relieve muscle soreness

Applications in the beauty field

- Slim legs
- Slimming

Circular Massage

Sequential compression massage from distal to proximal cyclically

Let everyone can own a suitable recovery device

• 12-chamber Air Pressure Therapy System

- Mode: A,B,C,D,E,F,G,H
- Pressure Range: 30-240mmHg
- Time: Optional treating time 5-99 minutes

► Machine Parameters

Voltage: AC220V or AC110V
Weight: 4.6 Kgs
Size: 30.4*25.4*16 cm

► Optional 5 Accessories





• Air Pressure Therapy System series

- Plug-in Air Pressure Therapy System (05~14)
- Rechargeable Air Pressure Therapy System (15~26)
- Varicose Vein / Deep Vein Thrombosis Therapy Device (27~28)

Let everyone can own a suitable recovery device

• 4-chamber Manual Version

- ⌚ Pressure Range: 30-240mmHg
- ⌚ Time: Optional treating time 0-30 minutes or run continuously

▼ Machine Parameters

Voltage: AC220V or AC110V Weight: 1.95 Kgs
 Size: 24*19.8*11.2 cm

▼ Optional 5 Accessories



VU-IPC02



• 4-chamber Digital Version

- ⌚ Mode: A,B,C
- ⌚ Pressure Range: 30-240mmHg
- ⌚ Time: Optional treating time 10,20,30 minutes

► Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2.1 Kgs
 Size: 24*19.8*11.2 cm

• 4-chamber Easy Version

-  Pressure Range: 30-240mmHg
-  Time: Optional treating time 0~30 minutes or run continuously




▼ Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2.2 Kgs
 Size: 30.1*23.7*12.7 cm

▼ Optional 5 Accessories



• 4-chamber Premium Version

-  Mode: A,B,C,D,E,F
-  Pressure Range: 30-240mmHg
-  Time: Optional treating time 1~99 minutes

► Machine Parameters



Voltage: AC220V or AC110V
 Weight: 2.3 Kgs
 Size: 30.1*23.7*12.7 cm



▲ VU-NIPC02 / VU-NIPC02H(with handle)

• 6-chamber Easy Version



-  Pressure Range: 30-240mmHg
-  Time: Optional treating time 0~30 minutes or run continuously

► Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2.2Kgs
 Size: 30.1*23.7*12.7 cm

▲ Optional 5 Accessories



Leg massage

Quickly eliminate fatigue and relieve muscle tension.

• 6-chamber Premium Version



VU-IPC04 / VU-IPC04H (with handle)

- Mode: A,B,C,D,E,F
- Pressure Range: 30-240mmHg
- Time: Optional treating time 1~99 minutes

► Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2.5 Kgs
 Size: 30.1*23.7*12.7 cm

Let everyone can own a suitable recovery device

• 8-chamber Premium Version



VU-IPC06 / VU-IPC06H (with handle)

- Mode: A,B,C,D,E,F
- Pressure Range: 30-240mmHg
- Time: Optional treating time 1~99 minutes

► Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2.65 Kgs
 Size: 30.1*23.7*12.7 cm

▲ Optional 5 Accessories



Rechargeable Air Pressure Therapy System series

The main purpose

- The device is intended for indoor and outdoor use by athletes, fitness enthusiasts and patients in the treatment of primary lymphedema, edema following trauma and sports injuries, post-immobilization edema, venous insufficiency, lymphedema.

Applications in the medical field

- Prevention of deep vein thrombosis (DVT), venous insufficiency
- Various primary, secondary and mixed edema such as lymphedema and postpartum edema
- Fracture, femoral head necrosis, soft tissue injury and other bone injuries
- Peripheral neuritis and diabetic foot caused by diabetes
- Rehabilitation patients after surgery
- Paralyzed
- Patients with varicose veins
- Long-time standing, sedentary or bedridden patients

Applications in sports field

- Relax and relieve fatigue after exercise
- Relieve muscle soreness

Applications in the field of beauty

- Slim legs
- Slimming

• **6-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D

Pressure Range: 30-240mmHg

Time: Optional treating time 5-99 minutes

▼ **Machine Parameters**

Li-ion Battery: 5000mAh/ 10.8V

Weight: 1.55 Kgs

Size: 22*13.82*8.77 cm

▼ **Optional 5 Accessories**



VU-IPC04R / VU-IPC06R

• **6-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D,E,F

Pressure Range: 30-240mmHg

Time: Optional treating time 5-99 minutes

► **Machine Parameters**

Li-ion Battery: 5000mAh/ 10.8V

Weight: 1.6 Kgs

Size: 22.8*14.2*9.5 cm



**Handheld;
Compact and
convenient**

Quickly eliminate fatigue and relieve muscle tension.



• **4-chamber Handheld Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D,E,F,G

Pressure Range: 30-180mmHg

Time: APP settings

► **Machine Parameters**

Li-ion Battery: 1600mAh/ 11.1V

Weight: 0.58 Kg

Size: 25.4*6.8*6 cm

▼ **Optional Accessories**



VU-IPCM4



4-chamber pants



4-chamber leg sleeve



• Rechargeable Air Pressure Therapy System

• 12-chamber Rechargeable Air Pressure Therapy System

- Mode: A,B,C,D,E,F,G,H
- Pressure Range: 30-240mmHg
- Time: Optional treating time 5-99 minutes

► Machine Parameters

Li-ion Battery: 9600mAh/ 10.8V
 Weight: 2.3 Kgs
 Size: 25.6*19*10 cm

► Optional Accessories

12-chamber leg sleeve

12-chamber arm sleeve

12-chamber pants



VU-IPC12B

• 8-chamber Rechargeable Air Pressure Therapy System



VU-IPCA8 / VU-IPCA4 / VU-IPCA6

- Mode: A,B,C,D,E,F
- Pressure Range: 30-240mmHg
- Time: Optional treating time 5-99 minutes

► Machine Parameters

Li-ion Battery: 5000mAh/ 10.8V
 Weight: 1.7 Kgs
 Size: 24.5*14.2*9 cm

▲ Optional Accessories

• **4-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D
 Pressure Range: 30-240mmHg
 Time: Optional treating time 5-90 minutes

Machine Parameters
 Li-ion Battery: 2500mAh/ 11.1V
 Weight: 1.3 Kgs
 Size: 22.4*15*8.4 cm

Optional Accessories

VU-IPCS4

4-chamber waist sleeve
 4-chamber shorts
 4-chamber leg sleeve
 4-chamber arm sleeve
 4-chamber pants

• **8-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D,E,F
 Pressure Range: 30-240mmHg
 Time: Optional treating time 5-90 minutes

Machine Parameters
 Li-ion Battery: 2500mAh/ 11.1V
 Weight: 1.5 Kgs
 Size: 22.4*15*8.4 cm

Optional Accessories

VU-IPCS8

8-chamber waist sleeve
 8-chamber shorts
 8-chamber leg sleeve
 8-chamber arm sleeve
 8-chamber pants

• **6-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D
 Pressure Range: 30-240mmHg
 Time: Optional treating time 5-90 minutes

Machine Parameters
 Li-ion Battery: 2500mAh/ 11.1V
 Weight: 1.4 Kgs
 Size: 22.4*15*8.4 cm

Optional Accessories

VU-IPCS6

• **6-chamber Rechargeable Air Pressure Therapy System**

Mode: A,B,C,D,E,F
 Pressure Range: 30-240mmHg
 Time: Optional treating time 5-99 minutes

Machine Parameters
 Li-ion Battery: 5000mAh/ 10.8V
 Weight: 1.6 Kgs
 Size: 24.5*15.2*9.2 cm

Optional Accessories

VU-IPCUG6



Rechargeable Air Pressure Therapy System

Let everyone can own a suitable recovery device

• 8-chamber Rechargeable Air Pressure Therapy System



Mode: A,B,C,D,E,F

Pressure Range: 30-240mmHg

Time: Optional treating time
5-99 minutes

► Machine Parameters




Li-ion Battery: 5000mAh/ 10.8V

Weight: 1.72 Kgs

Size: 24.5*15.2*9.2 cm

▲ Optional Accessories

• **5-chamber Rechargeable Air Pressure Therapy System**

-  Mode: A,B,C,D,E,F
-  Pressure Range: 30-240mmHg
-  Time: Optional treating time 5-99 minutes

► **Machine Parameters**

Li-ion Battery: 2500mAh/ 11.1V
 Weight: 1.42 Kgs
 Size: 13.4*15*9.4 cm




► **Optional Accessories**

5-chamber leg sleeve
 5-chamber legs double tube



▲ VU-IPC5-V2 / VU-IPC5-V1

• **6-chamber Rechargeable Air Pressure Therapy System**

-  Mode: A,B,C,D,E,F
-  Pressure Range: 30-240mmHg
-  Time: Optional treating time 5-99 minutes

► **Machine Parameters**

Li-ion Battery: 2500mAh/ 11.1V
 Weight: 1.5 Kgs
 Size: 13.4*15*9.4 cm

**Handle design,
 convenient for carrying**
 Circulating massage, relaxes the body



▼ **Optional Accessories**



▲ VU-IPC6-V2 / VU-IPC6-V1







Arm and leg massage

Kneading, pressing,
relieving physical and mental fatigue

• Varicose Vein / Deep Vein Thrombosis Therapy Device



VU-IPC05C

-  Number of modes: 1 modes
-  Number of chambers: 2 chambers
-  Pressure Range: 0-60mmHg
-  Time: continuous work

► Machine Parameters

Li-ion Battery: 3000mAh/ 3.7V
Weight: 0.2088 Kg
Size: 14.3*7.1*2.7 cm

► Standard Unit







sleeve



charger

• Varicose Vein / Deep Vein Thrombosis Therapy Device

-  Number of modes: 2 modes
-  Number of chambers: 2 chambers
-  Pressure Range: 0-90mmHg
-  Time: continuous work

► Machine Parameters

Li-ion Battery: 3000mAh/ 3.7V
Weight: 0.2088 Kg
Size: 14.3*7.1*2.7 cm



VU-IPC05B / VU-IPC05B-V2

Intended use

· This equipment is intended for use by healthcare professionals to use on patients or for patients to use at home by themselves, this device by helping to stimulate blood flow in the legs, help prevent varicose veins, as well as deep vein thrombosis (DVT) formation.

Medical Use

- Help prevent deep vein thrombosis (DVT); improve blood circulation; reduce postoperative pain and swelling; shorten wound healing time;
- Help treat and heal stasis dermatitis, venous eczema, stasis ulcers, arterial and diabetic leg ulcers, chronic venous insufficiency, and edema reduction;
- It can also be used to prevent deep vein thrombosis in people who sit or lie down for long periods of time.



• Hand Rehabilitation Training Device

• Hand Rehabilitation Training Device

► Machine Parameters

Mode: A/B/C/D
 Time: Optional treating time 5-60 minutes
 Main interface: Type-C charging port
 Li-ion Battery: 6000mAh/ 7.4V
 Size: 20.6*14.6*8.5cm

► Accessories List



Recovery glove Mirror glove Air glove Boot Recovery ball Charger



▲ VU-FR01

• Hand Rehabilitation Training Device

► Machine Parameters

Time: Optional treating time 0-30 minutes
 or continuous work
 Voltage: 110V~220V
 Size: 30.1*23.7*12.7cm

► Accessories List



Recovery glove Mirror glove Recovery ball



▲ VU-FR02

• Hand Rehabilitation Training Device

► Machine Parameters

Time: Optional treating time 0-30 minutes
 or continuous work
 Li-ion Battery: 3200mAh/ 7.4V
 Size: 16.2*8.56*4.38 cm

► Accessories List



Recovery glove Mirror glove recovery ball



▲ VU-FR03



• Anti-Bedsores Inflatable Mattress With Pump

• Anti-Bedsores Inflatable Mattress With Pump

- Mode: turn left, turn right, back up, lie down
- Number of chambers: 5 chambers
- Pressure Range: 20-120mmHg
- Time: continuous work

► Machine Parameters

Voltage: AC220V or AC110V
 Weight : 1.4 Kgs
 Size : 24*19.8*11.2 cm

► Accessory



▲ VU-AMP01

• Anti-Bedsores Inflatable Mattress With Pump

- Mode: turn left, turn right, turn left&right, sit up, static/wave, (leg bend)
- Chambers: 5 (Without leg bend)
6 (With leg bend)
- Pressure Range: 20-110mmHg
- Time: continuous work

► Machine Parameters

Voltage: AC220V or AC110V
 Weight: 2 Kgs
 Size: 30.1*23.7*12.7 cm



▲ VU-AMP02



• Medium-frequency Massager

Let everyone can own a suitable recovery device

• Medium-frequency Massager

► Machine Parameters

Mode: kneading massage/pounding massage/
acupressure massage

Time: Optional treating time 5-30 minutes

Main interface: Type-C charging port

Li-ion Battery: 230mAh/ 3.7V

Power: 4 VA

Size: 5*5*1.1 cm

► Accessories List

Massage pad, Charging wire



▲
VU-MF01



Main purpose

Muscle fatigue such as shoulder weakness, fatigue after exercise, unilateral or bilateral shoulder discomfort, neck and shoulder fatigue, waist fatigue, etc.



• Migraine Headache Relief

• Migraine Headache Relief



Host and
Massage pad



Bluetooth



continuous work

► Machine Parameters

Mode: relaxing/headache massage/sleeping

Time: Optional treating time 5-30 minutes

Main interface: Type-C charging port

Li-ion Battery: 80mAh/ 3.7V

Power: 0.37 VA

Size: 3.7*3.7*1.2 cm

► Accessories List

Massage pad、Charging wire



▲
VU- HR01

Main purpose

Treatment of headaches, migraines, tension induced insomnia, and auxiliary treatment of anxiety and depression.



• Medical Parallel Transport Bed

• Structure composition

Bed Frame, adjustable laying board on wheels, lifting mechanism, sliding and stretching mechanism, wheels, guardrail, control device, handheld controller, brake, handlebar, battery.



• Main functions

- a) Bed surface raise or down
- b) Bed surface extension board can move horizontally
- c) Bed surface transmission cloth can rotate left and right
- d) Guardrails can be raised or down
- e) Equipped with a foot wheel locking device
- f) Equipped with an electric power-assisted wheel



• Basic Parameters

Operating mode	-----	Intermittent operation
Grounding condition	-----	Grounded
Input voltage	-----	AC100-240V, 50Hz
Input power	-----	250 VA
Dimensions	-----	220*70*68 cm





• Electric Medical Bed/nursing Bed

• Product advantages

- 1、 Built-in clinical nursing features, optimize nursing process, improve nursing efficiency, and reduce common risks.
- 2、 The unique three-fold turning function improves the quality of life.
- 3、 The configuration of the human-computer interaction system, the online display of the clinical characteristics and real-time status of the patient or the elderly.

• Features

- 1、 Initial state
- 2、 Automatic left/right turning function
- 3、 Automatic back rest, knee rest functions, and one-button molding

An electric multifunctional medical/nursing bed designed to provide a safe environment for the care and rehabilitation of patients and the elderly.





• Rolling Massager

► Machine Parameters

Input Voltage: DC 24V
Power: 70 w
Weight: 4.5 kgs
Size: 48*24*18.5 cm

► Main Features

1. Ten rollers rotate and roll massage
2. Three-level strength design (users can freely adjust the appropriate roller speed)
3. It can massage the soles of the feet, calves, and waist flexibly

• Rolling Massager

Intended use

Rolling massager is to use rollers to repeatedly roll on the surface of the feet and waist, promote local blood circulation, enhance body metabolism, and achieve the goal of relaxing the body and mind.





