

Hubei YJT Technology Co., Ltd



Hubei YJT Technology Co., Ltd

Email: weldon@yjt.com.cn

Tel: +86-27-87646661

Website: www.yjthnc.com
www.yjtmedical.com

Add: Floor 8, Block 7, No.62 Optical Valley Ave, East Lake High-tech
Development Zone, Wuhan, China



Help your health, light your life



Real Review, Real Partner

It was very nice deal with this company. A good product locking for more deals with them in future.

———— **Natalia Burtman**

YJT is a reliable company. It always keep their product in good quality level. We make business with them over 6 six years and never make us disappointed. They also provide good after sales services and always have a prompt response. We are satisfied with them in services and products.

———— **Frank Johnson**

I am very happy to have gone through a wonderful five years with YJT. Your professionalism in products and excellence in quality are the cornerstones of our cooperation. It is a great honor to witness the 18th birthday of YJT. I believe we can create more and more innovative products to serve the society in future.

———— **Andrew Smith**



Company Profile

- Hubei YJT Technology Co., Ltd. was established in 2005.
- Headquarter located in Wuhan China with 3,000 square meters office.
- Industrial park located in Ezhou covered more than 57,000 square meters.
- With 18 years experience and deep understanding of physiotherapy, we are the first listed company (stock number 870130) in the phototherapy industry.
- Own up to 50 qualified engineers and 10 automatic producing line.
- Products served to 1.5 million families around the world.
- Expertise in physiotherapy industry
 - Research and development of innovative product.
 - Manufacturing and marketing of physical therapy device.
 - OEM and ODM service
- Successfully performed 800 manufacturing service to US, Europe, Asia etc.
- Approved by FDA, CE, ISO9001, ISO13485, RoHs.



Our Statement



Mission

Keep your life away from pain

Vision

Strive for Health of Millions of Families

Slogan

Help your health, light your life



President's Oration



YJT will lead the way to the future of phototherapy industry, not only engage ourself in supplying qualified products and high class service, we also would like every single device from us could help people around the world to relieve pain, keep away from high blood pressure, solve the sub-health problem, offer efficient treatment solution etc.

Through the commitment to our products quality and vision, we will make constant innovation to exceed our clients' expectation. We expect people still remember YJT after decades later even we want to be an enterprise of hundred years.


Let us strive for health of millions of families.

President of YJT



Certificates

FDA registered factory, approved by ISO9001 and ISO13485, certified by CE, 510K, RoHs.



Main Product Range



Handy laser therapy device for personal pain relief



Nail fungus laser therapy device for onychomycosis

Wrist type laser therapy device for high blood viscosity control



Laser comb and laser cap for hair regrowth




Class IV laser therapy device for pain relief in clinic



LED light therapy device for male and female private problem

Millimeter wave therapy device for diabetes treatment



High potential therapy device for sub-health problem



Handy laser therapy device



Item	HD-Cure	HD-CureS	HD-Cure Pro	HD-CurePro 01
Laser wavelength	650nm&808nm	650nm&808nm	650nm&808nm	808nm&650nm
Laser diodes	12*650nm,3*808nm	12*650nm,4*808nm	6*650nm,6*808nm	5*808nm ,12*650nm
TENS power	×	×	40-200Hz	×
Laser power	510mW	660mW	930mW	810mW
Laser working mode	Continues and Pulse	Continues	Continues	Continues and Pulse
Battery	5200mAh	2600mAh	2600mAh	3350mAh
Comb cover	×	×	×	✓

Indication:

Pain relief
Knee arthritis
Joints and bones injuries
Sports injuries & sprain
Soft tissue injuries
Acupuncture rehabilitation therapy
Athletics system injuries
Wound&Ulcer acupuncture
Anti-inflammation
Hair growth(HD-CurePro 01)

Advantage:

- FDA registration
- Suitable for personal health care
- Physiotherapy device, no side effect, no cross infection, non- invasive, natural healing
- Portable design,easy to operate and use at anywhere anytime
- Large lithium battery capacity, long battery life



Multi-functional laser therapy machine



Item	Description
Laser medium	GaAlAs Semiconductor
Laser wavelength	808nm and 650 nm
Terminal laser output	Probe A:1*808nm Probe B:4* 808nm and 11* 650nm
Terminal maximum power output	Probe A is 180mW and Probe B is 775mW
Laser power level	5 levels

Indication:

Pain relief of acute/chronic Injure
Wound healing
Fracture healing
Anti-inflammation and rehabilitation therapy

Advantage:

- Perfect combination of medical laser 808nm and home laser 650nm
- Security and stability of home laser as well as the exact effect of medical laser
- Available from 5mw to 775mw and can be used while still plugged
- Creative systemic treatment mode of large surface exposure and acupuncture point exposure to irradiation
- A laser stabilizer specially added to perform with more stability



Class IV laser therapy device



Item	Description
Laser medium	GaAlAs Semiconductor
Laser wavelength	808nm and 650nm
Terminal laser output	4 × 808nm laser beams, 500mW each beam
	10 × 650nm laser beams, 30mW each beam
Total output power	2300mW
Battery capacity	Built-in lithium battery 5200mAh
Working mode	Continuous and pulsed mode 10Hz
Case Material	Aluminum
Display	White LED display
Time setting	30s-3 minutes options, each grade increases 30 seconds

Indication:

Pain relief of Acute/chronic injury
Wound healing
Fracture healing
Anti-inflammation and rehabilitation therapy

Advantage:

- Class IV laser
- Perfect Combination of 650nm laser, 808nm infrared laser
- Physiotherapy device, no side effect, no cross infection, non-invasive, natural healing
- Good therapy effect, clinical trail proved
- Key switch and emergency button for double security



Nail fungus laser therapy device



Item	HZJ-01	HZJ-02	HZJ-03
Laser wavelength	905nm	905nm	980nm
Blue light wavelength	465nm	465nm	465nm
UV wavelength	×	×	270nm
Battery	1200mAh	1200mAh	900mAh
Working time	7 minute	7 minute	7 minute
Digital display	N/A	Yes	No

Indication:

Nail/Toe nail fungus infections
Onychomycosis

Advantage:

- First manufacturer in China
- Non-invasive, no side effects, painless, pure physiotherapy
- Portable design, rechargeable lithium battery inside
- Easy for using at home and in a travelling



Wrist type laser therapy device



Item number	SL-03	SL-04	SL-08
Laser wavelength	650nm	650nm&808nm	650nm
Laser diodes	18*650nm	12*650nm,1*808nm 2*465nm blue LED 2*586nm yellow LED 2*520nm green LED	12*650nm
Heart rate	×	×	yes
Laser power	90mW	210mW	60mW
Laser working mode	Continues and Pulse	Continues and Pulse	Continues and Pulse
Battery	1900mAh	1580mAh	1580mAh

Indication:

Hyperviscosity
Hyperlipemia
Hypertension
Diabetes
Rhinitis
Ischemic cardiac-cerebral vascular diseases

Advantage:

- High quality imported laser head, long life-span
- Humanistic and scientific design to fit human physiological curve and ergonomics
- Medical use high capacity lithium battery for the energy -saving
- High-class and durable metal enclosure which is conducive to heat diffusion for the power stable
- Two output modes: pulse and continues, which have manual and automatic switch available



Laser comb for hair regrowth



Model number	Lasercomb-001	Lasercomb-002
Number of light guides	7pcs	9pcs
Laser Wavelength	650nm	650nm
Adapter voltage	DC5V	DC5V
Battery	2600mAh	2600mAh
Massaging	Vibration	Vibration

Indication:

Anti-hair loss
Hair regrowth
Scalp massaging

Advantage:

- Laser and massage all in one machine, combined light irradiation with massage mode
- With 650nm red laser light
- Light guide teeth, intensive layout, more focused illumination surface
- Mini size, portable for home use anytime anywhere
- Power switch and massage key in elegant design



Laser cap for hair growth



Model		Hat-01
Laser Medium	GaAIAS Semiconductor	GaAIAS Semiconductor
Number of laser diodes	118pcs laser	82pcs laser
Wavelength	650nm	650nm
Voltage	100-240V	100-240V
Battery	No	1900mAh x 2pcs
Working Mode	Continuous or pulse mode	Continuous or pulse mode
Time	15mins or 30mins adjustable	30mins

Indication:

Treatment of hair loss
Hair regrowth and repair damaged hair
Balance scalp oil secretion
Improve the quality of the hair

Advantage:

- Low energy laser, revitalize and thicken hair
- laser diodes---fully covering the entire scalp
- Free hands treatment, use worry-free in multiple scenarios
- Directly reach the subcutaneous tissue 6-8mm, relax the scalp
- Safe to use, no side-effect and ergonomics design



Male prostate red and blue light therapy machine



Item	Description
Red chips	630nm \pm 10nm
Blue chips	465nm \pm 10nm
Power intensity each chip	10mW/cm ²
Timing Function	5-60 minutes adjusted Blue light time suggested 5 minutes Red light time suggested 30 minutes
Rated power input	500VA

Indication:

Benign prostate hyperplasia
Perineal pain
Difficulty urinating
Frequent urination
Nocturnal polyuria

Advantage:

- Intelligent temperature control chip, accurate and safe
- HD Touch Screen
- Blue light is sterilization, red light is anti-inflammatory
- Wearable structure, easy to use
- Free replacement for after sales service



Gynaecology red and blue light therapy device



Item	Description
Red light wavelength	630nm ± 10nm
Blue light wavelength	465nm ± 10nm
Power intensity each diode	10mW/cm ²
Timing Function	0- 60minutes adjustable Error is controlled within 10% A bee sound means the treatment and red light & ozone output stop
Light power level	3 levels

Indication:

Mild cervical erosion
Mild non-specific vaginitis
Cervicitis
Pelvic inflamma-tory disease

Advantage:

- Blue light can kill bacteria, red light can reduce inflammation
- With extra machine host to operate the led probes, more upscale for clinic and health centre treatment
- Free replacement for after sales service



Knee massager with heating



Thermal therapy	LED-K
Number of light guides	70pcs (LED infrared light diodes)
LED Light Wavelength	630nm,660nm,850nm
Thermal therapy	heating,2 modes
Massaging	Vibration, 2 modes
Adapter voltage	DC5V
Battery	2450mAh

Indication:

Wounds healing
Pain relief
Anti-inflammation
Fibromyalgia
Nerve injuries

Sports injuries
Muscle injury
Arthritis
Massager

Advantage:

- Design for joint pain,wearable
- 3 in 1 multi-functional
- No side effects
- Easy operation at home
- Safe for all skin types and ages



Millimeter wave therapy device



Model number	MMW-01	MMW-02
Wavelength	Millimeter-wave 7.5mm-10mm	
Output power density	4mW/cm ² ~ 9mW/cm ²	
Operating frequency	30GHz ~ 45GHz	
Rated input power	<30VA	<35VA
Probes	1 probe	2 probes
Screen	High definition and press panel screen	High definition and touch panel screen
Strap	✓	✓

Indication:

Diabetic and complication diseases, such as diabetic foot
Diabetic nephrosis, etc. Inhibiting cancer&tumor cells growth, and relieve the pain

Advantage:

- Non-thermal effect: no thermal, no pain, no stimulation, physical therapy
- The cumulative effect of time: with the increase in the time of treatments, the treatment effect is better
- High frequency electromagnetic wave with 1-10mm wavelength
- More than 95% efficacy in diabetic complications



Electrical potential therapy equipment



Model number	HPT	HPT001	DW9202-02	HPT003
High potential output voltage	5000V 7000V 9000V		8000V 10000V 12000V	9000V 12000V 14000V
Negative potential output voltage	-4000V -5000V -6000V	-4000V -5500V -7000V	-4000V -5500V -7000V	-6000V -8000V -10000V
If potential	80Khz			
Intelligent mode	5mins gap intelligence circulation			
Output frequency	50Hz/60Hz			
Pain releasing handle	✓	✓	✓	✓
Negative ions	×	×	✓	×
Warm moxibustion	×	✓	×	✓

Indication:

Headache, insomnia, bone&joint pain, chronic constipation

Advantage:

- By stimulating and reconstructing the natural magnetic field through potential, the treated people can return to nature and balance the various functions of the body
 - Repair various disorders of the human body caused by the destruction of the natural environment
 - Supplement the static resting energy of the human body and eliminate free radicals
 - Safe and stable, with overheat protection and output voltage stabilization device
- Multiple modes of treatment in one



CES insomnia therapy device



Model number	HSM-01	HSM-02
Pulse modulation frequency range	15HZ	
Battery Capacity:	1200mAh	800mAh
charging time:	4-5 hours	
Ear clip	2 clips	
TENS	2 electrode pads	
Screen	LCD	LED
Cotton flakes	✓	✓

Indication:

Anti-insomnia, asleep faster, anti-depressant, anti-anxiety, relax&calm, pain relief with TENS pad

Advantage:

- No pain, no side effects
- Equipped with TENS body acupoint electrode to relieve pain
- Rechargeable lithium battery inside
- Adjustable treating time and power
- Working in low frequency pulse



Integrated treatment instrument for myopia and amblyopia



Item	Description
Red light wavelength	650nm ± 20nm
Number of laser diodes	Eye laser utilize 2pcs + Acupoint laser utilize 20pcs
Running mode	Continuous
Internal power supply	DC3.7V
Input	100V-240V;50/60HZ;0.5A-0.2A
Laser output power	2mW ± 25%

Indication:

Amblyopia
 Refractive and inflammatory amblyopia,
 Refractive amblyopia,
 Strabismus amblyopia
 Myopia patients, are only used for people aged 7-18

Advantage:

- Combine laser with acupoint treatment
- Portable & Rechargeable
- Non-invasive, no side effects, painless, pure physiotherapy



Air compression leg massager



Item	Description
Dimensions	8.7" × 2.75" × 2" (221mm × 70mm × 48.3mm)
Weight	Approx × 0.81lb (368g)
Battery capacity	Built-in lithium battery 2600mAh
Working mode	Mode1, Mode2, Mode3, Mode4
Source of power	90-264VAC, 50-60HZ
Mode of operation	Continuous

Indication:

Aid in the prevention of DVT
 Enhance blood circulation
 Diminish post-operative pain and swelling
 Reduce wound healing time
 Aid in the treatment and healing of :stasis dermatitis, venous stasis ulcers, arterial and diabetic leg ulcers, chronic venous insufficiency and reduction of edema in the lower limbs.

Advantage:

- Multiple modes can be switched
- Portable & Rechargeable
- Non-invasive, no side effects, painless, pure physiotherapy



Nano belt for waist pain relief



Item	Description
Dimensions	1450mm × 240mm × 40mm
Vibration pattern	3pattern
Heating level	3levels
Battery capacity	2600mAh
Using time	20 min a time
Source of power	110V-240V, 50-60HZ

Indication:

Waist pain relief
 Waist massage
 Muscle injury
 Back supporting
 Improve blood flow

Advantage:

- Wireless & Portable & Rechargeable
- Non-invasive, no side effects, painless, pure physiotherapy



ODM Service



